



Sample cells with caps for 2100 portable turbidimeters, 10 mL, 6 pcs

Nr. produs:	2434706
RON Preț (fără TVA):	Contact
Disponibil	

Replacement sample cells for use with portable turbidimeters. Pack of six round sample cells with caps.

Specificații

Cantitate:	6 /pk
Greutate:	0,1 kg
Material:	Glass
Nota de subsol:	For use with Hach 2100 Series Portable Turbidimeters
Scop alimentare:	6 sample cells with caps
Volum probă:	10 mL