



Be Right™



## pH test kit, 6.5 - 8.5 pH, model 17H

Nr. produs: 147008

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Simple, low-cost color disc method is the most accurate visual test kit method.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specificații

Domeniu de măsurare:	6.6 - 8.4 pH
Domeniu de masurare :	NA
Greutate bruta:	1
Model:	17H
Model test kit:	17H
Numar teste:	200
Nume metoda:	Color disc/Phenol red
Parametru:	pH, mid range
Pauza:	NA
Platforma:	Disc culoare
Smallest increment:	0.1 pH
Tipul carcasei:	B