



Iron color disc test kit, model IR-18B

Nr. produs:

146401

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	High Single Range
Interval de măsurare standard:	0.25 - 7 mg/L Fe
Model:	IR-18B
Model test kit:	IR-18B
Numar teste:	100
Nume metoda:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron
Parametru:	Iron, medium range - as Fe
Platforma:	Disc culoare
Smallest increment:	0.25
Tip:	Simplu Metric
Tipul carcasei:	B