



Iron color disc test kit, model IR-18A

Nr. produs:

146500

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	Low Single Range
Interval de măsurare standard:	0.02 - 1.00 mg/L Fe
Model:	IR-18A
Model test kit:	IR-18A
Numar teste:	100
Nume metoda:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron
Parametru:	Iron, medium range - As Fe**
Platforma:	Disc culoare
Smallest increment:	0.02
Tip:	Simplu Metric
Tipul carcasei:	B