



## Testkit, glycol visual, model EG-1

**Nr. produs:** 2186400  
**RON Preț (fără TVA):** Contact  
**Nu este afișată data expirării**

Produse periculoase



Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specificații

Domeniu de măsurare:	Presence or absence of glycol, 0- greater than 150 mg/L
Model test kit:	EG-1
Numar teste:	25
Nume metoda:	Visual/  Appearance of color
Parametru:	Glycol - Ethylene glycol  in oil or water
Tipul carcasei:	G