



Test kit, manganese color disc, model MN-PAN

Nr. produs:

2350800

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	0 - 0.70 mg/L Mn
Interval de măsurare standard:	0 - 0.70 mg/L Mn
Model:	MN-PAN
Model test kit:	MN-PAN
Numar teste:	50
Nume metoda:	Color disc/PAN
Parameter:	Manganese
Parametru:	Manganese - As Mn
Platforma:	Disc culoare
Smallest increment:	0.05
Tipul carcasei:	D