



## Test kit, total chlorine, model CN-DT, 20-2000 mg/L

**Nr. produs:** 2471100

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Produse periculoase



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Kit contains all required reagents and apparatus to perform analysis in the field**

### Specificații

Domeniu de măsurare:	20 - 2,000 mg/L $\text{Cl}_2$
Model test kit:	CN-DT
Nota de subsol:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$ , ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Numar teste:	100
Nume metoda:	Digital Titrator/Thiosulfate
Parametru:	Chlorine, mid and high range - As Total $\text{Cl}_2$
Platforma:	Titrator digital
Smallest increment:	0.2 - 5.0
Tipul carcasei:	F