



Sulfite Test Kit, Model SU-5

Nr. produs: 148002
RON Preț (fără TVA): Contact
Disponibil

Produse periculoase



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	1.0 - 20 mg/L & 10 - 200 mg/L SO_3^{2-}
Model test kit:	SU-5
Numar teste:	100
Nume metoda:	Drop count, titration Iodometric
Parametru:	Sulfite
Platforma:	Drop Count Titration
Smallest increment:	1 10
Tipul carcasei:	D