



## Test kit, alkalinity, model AL-AP

**Nr. produs:** 2444300  
**RON Preț (fără TVA):** Contact  
**Expediere în 2 săptămâni**

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specificații

Domeniu de măsurare:	- 8 gpg
Domeniu de masurare :	1 - 20 gpg**
Greutate bruta:	1
Model test kit:	AL-AP
Nota de subsol:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$ , ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.  **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Numar teste:	100
Nume metoda:	Drop count titration/Sulfuric acid
Parametru:	Alkalinity - (P) & Total (MO) as $\text{CaCO}_3$
Platforma:	Titration
Smallest increment:	0.4 gpg  1 gpg
Tipul carcusei:	B