



Test kit, alkalinity, model AL-AP

Nr. produs:

2444300

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	0 - 8 gpg & 1 - 20 gpg
Model test kit:	AL-AP
Nota de subsol:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Numar teste:	100
Nume metoda:	Drop count titration/Sulfuric acid
Parametru:	Alkalinity - (P) & Total (MO) As CaCO ₃
Platforma:	Drop Count Titration
Smallest increment:	0.4 gpg
	1 gpg
Tipul carcasei:	B