



Test kit, ozone, model OZ-2

Nr. produs: 2064400
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Produse periculoase



Hach's most economical test kit for ozone measurement utilizes the simple, low-cost color disc method.

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 2.2 mg/L O ₃
Interval de măsurare standard:	0 - 2.2 mg/L O ₃
Model:	OZ-2
Model test kit:	OZ-2
Numar teste:	100
Nume metoda:	Color disc/DPD
Parameter:	Ozone
Parametru:	Ozone - As O ₃
Platforma:	Disc culoare
Smallest increment:	0.05
Tipul carcasei:	B