



Bottle, BOD, 300 mL, 6/pk

Nr. produs: 62106
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Set of 6 glass BOD bottles, 300 mL, with stoppers. For use with LBOD sensor to determine the Biological Oxygen Demand. Borosilicate clear glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included.

Specificații

Body Material:	Glass
Cantitate:	6 /pk
Capacitate:	300 mL
Material:	Borosilicate Glass
Numar metoda:	8043
Package size UOM:	6pk