



Bottle, BOD, 300 mL, 24/pk

Nr. produs: 62124
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

Specificații

Body Material:	Glass
Cantitate:	24 /pk
Capacitate:	300 mL
Material:	Borosilicate Glass
Package size UOM:	24pk