



## Beaker, polypropylene, low form, 400 mL, 6/pk

Nr. produs: 108078

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specificații

Capacitate:	400 mL
Diametru exterior :	99.1 mm with spout
Material:	Polypropylene, Low Form
Material corp:	Polypropylene
Unitate:	6 /pk