



Be Right™



## Cylinder, graduated, polypropylene, 50 mL $\pm 0.6$ mL, 1.0 mL divisions

Nr. produs: 108141

RON Preț (fără TVA): Contact

Disponibil

Hach offers graduated cylinders in a variety of classes and sizes for liquid transfer or mixing. Choose Class A graduated cylinders for the best tolerance limits.

Volumetric glass cylinders are listed as "To Contain" or "To Deliver". Glassware designed to contain are typically marked with TC. When liquid is poured out of the TC cylinders, a small amount remains behind, clinging to the sides of the glass. It will deliver slightly less than what was measured. To deliver cylinders are typically marked TD. When liquid is poured out of the TD cylinders, a small amount of liquid is left clinging to the walls, but the TD design accounts for this and the full measured amount is what is delivered. Pipets and burets are also designed to deliver.

---

### Specificații

Cantitate:	1 each
Capacitate:	50 mL
Certificări:	None
Material corp:	Polypropylene
Metoda de masurare:	To Deliver/To Contain
Tip:	Simple Metric