



CEL/890 Advanced portable laboratory

Nr. produs:

2994000

Produse periculoase

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.



Offering on-site results in an easy-to-use, rugged, portable lab, Hach's Oil & Gas Water Quality Laboratory covers parameters critical to effective source water and base water analysis. Depending on water composition, this kit may also be used to support analysis of produced water, flowback water, water treatment, and enhanced oil recovery applications.

Specificații

Aplicatie:

Fracing

Conținutul cutiei:

This laboratory includes: DR 900 Colorimeter, HQ40d Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 critical parameters, apparatus and accessories, and two convenient, rugged carrying cases.

Conținutul cutiei

This laboratory includes: DR 900 Colorimeter, HQ40d Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 critical parameters, apparatus and accessories, and two convenient, rugged carrying cases.