



Be Right™



Rectangular sample cell, 1"x1", 25 mL, glass (2 pieces)

Nr. produs: 2095000

RON Preț (fără TVA): Contact

Disponibil

Accessories for spectrophotometers/spectral measurements of liquids.

For use with the DR 2010 and DR 2000 spectrophotometers. Comes as a set of matched pair cells.

Accessories for spectrophotometers/spectral measurements of liquids.

For use with the DR 2010 and DR 2000 spectrophotometers. Comes as a set of matched pair cells.

For use with the DR/2010 and DR/2000 Spectrophotometers. Comes as a set of matched pair cells.

Specificații

Cantitate:	2
Lungime de unda (optica):	1 inch rectangular, 25 mL
Material:	Glass