



DR5000 Multi-cell holder

Nr. produs: A23618

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Accessories for spectrophotometers/spectral measurements of liquids. The holder is designed for 10, 20 or 50mm rectangular cells, or 1-inch round and rectangular ones.

Accessories for spectrophotometers/spectral measurements of liquids. The holder is designed for 10, 20 or 50mm rectangular cells, or 1-inch round and rectangular ones.

The DR 5000 Spectrophotometer comes equipped with an adapter that can accommodate the following cuvette types: ? 10 / 20 / 50 mm rectangular cells ? 1-in. round cell / AccuVac« ? 1-in. square cell

Specificații

Cantitate:	1 pcs
Dimensiuni (Înălțime x Lățime x Adâncime):	42 mm x 107 mm x 107 mm
Greutate:	0,094 kg
Material:	ABS