



Be Right™



## POCKET Colorimeter II Sistem Pocket Colorimeter™ II - OD

Nr. produs: 5870003

Produse periculoase

**ARTICOL ÎNVECHIT**  
Acest articol nu mai este disponibil.



Portable Colorimeter programmed for the determination of Dissolved Oxygen, HRDO method. In carrying case, complete with reagents, cuvettes and operating instructions.

Portable Colorimeter programmed for the determination of Dissolved Oxygen, HRDO method. In carrying case, complete with reagents, cuvettes and operating instructions.

### Specificații

Acuratețe lungime de und#:	Lungime de undă fixă $\pm 2$ nm variază în funcție de model
Carcasă impermeabilă (IP) :	IP67
Certificarea conformității:	CE
Detector:	Silicon detector
Dimensiuni (Înălțime x Lățime x Adâncime):	155 mm x 61 mm x 35 mm
Domeniu de măsurare:	0.2 - 10.0 mg/L O <sub>2</sub>
Garanție:	2 ani
Greutate:	0,23 kg
Instrument:	Kitul include instrumentul PC II
Lamp Type:	LED (diodă emițătoare de lumină)
Limbi manual:	...
Lungime de unda:	Fixed wavelength $\pm 2$ nm varies with model
Monitor:	LCD, cu iluminare de fundal
Numar teste:	25
Nume metoda:	Colorimetru/  AccuVac®
Parametru:	Oxigen dizolvat (OD)
Smallest increment:	0.1

### Conținutul cutiei

Sistem de analiză Pocket Colorimeter II, oxigen dizolvat - între 0,2 și 10,0 mg/l, cu set de reactivi (25 de teste). Toate instrumentele Pocket Colorimeter II sunt furnizate în cadrul unui kit gata de utilizare, care include un set de reactivi (doze unități de reactivi premăsurate), cutie de transport solidă personalizată și manual.