



pH electrode 5303, gel electrolyte, 3x diaphragms, 0-80 °C, 6 bar, S8 screw head

Nr. produs: LZU5303.97.0002

RON Preț (fără TVA): Contact
Disponibil

pH electrode with gel electrolyte and S8 screw head. 0 - 14 pH, 0 - 80 °C, up to 6 bar. 3x ceramic diaphragms.

Applications: swimming pools, drinking water, gas

scrubbing towers, boiler water, detoxification,

wwtp outlet water, etc.

Specificații

Conexiunea cablului:	S8
Diafragma:	3 ceramic
Diametru:	12 mm
Domeniu de temperatură:	0 - 80 °C
Electrolit:	Gel
Garanție:	6 luni
Interval de presiune:	până la 6 bar
Lungime:	120 mm
Material:	Sensor: glass
Measuring range pH:	0 - 14 pH
Parametru:	pH
Reference system:	Encapsulated Ag/AgCl crystals