



Molybdenum Reagent Set, Powder Pillows, 0-35.0 mg/L Mo

Nr. produs:

2243400

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Powder reagents in sealed packages for the determination of Molybdenum, Colorimetric Mercaptoacetic Acid Method. For 25ml sample.

Powder reagents in sealed packages for the determination of Molybdenum, Colorimetric Mercaptoacetic Acid Method. For 25ml sample.

For determination of Molybdenum by the Mercaptoacetic Acid Method (Hach method 8036). Range: Up to 35.0 mg/L as Mo⁶⁺. Sample Size: 25 mL. Set Contains: MolyVer[®] 1, 2 and 3 Molybdenum Reagent Powder Pillows for 100 tests.

Specificații

Domeniu de măsurare:	0 - 35,0 mg/L Mo
Domeniu de masurare :	25 mL
Instrument:	DR/3000, DR/2000
Numar metoda:	8036
Numar teste:	100
Nume metoda:	Mercaptoacetic Acid
Platforma:	Pernu#e cu pudră