



Reagent powder pillows, molybdenum 1, 50 mL, pk/100

Nr. produs: 2352769

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Powder reagent in sealed package for the determination of Molybdenum, 0.02-3.00mg/l Mo.

Powder reagent in sealed package for the determination of Molybdenum, 0.02-3.00mg/l Mo.

For determination of low range molybdenum by the ternary complex method (Hach method 8169). Range: 0.02-3.00 mg/L Mo Sample Size: 50 mL

Specificații

Cantitate:	100 /pk
Numar teste:	100
Nume metoda:	Ternary Complex
Platforma:	Pernu#e cu pudră