



Porphyrin 2 Reagent Powder Pillows, 25 mL, pk/100

Nr. produs:

2187569

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



Powder reagent in sealed package for the determination of Copper, Bicinchoninate Method, 2-210 µg/l Cu. For 25ml sample.

Powder reagent in sealed package for the determination of Copper, Bicinchoninate Method, 2-210 µg/l Cu. For 25ml sample.

For determination of copper by the Porphyrin method. For 25 mL samples

Specificații

Cantitate:	100 /pk
Nume metoda:	Porphyrin
Platforma:	Pernu#e cu pudră