



Reagent set, copper, porphyrin, 25 mL

Nr. produs:

2242700

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Powder reagents in sealed packages for the determination of Copper, Porphyrin Method. For 25ml sample.

Powder reagents in sealed packages for the determination of Copper, Porphyrin Method. For 25ml sample.

For determination of copper by the Porphyrin Method. Hach Method 8143. Range: 2 - 210 µg/L Sample Size: 25 mL. Reagent Set includes: Copper Masking Reagent Powder Pillows, Porphyrin 1 and 2 Reagent Powder Pillows for 100 tests.

Specificații

Domeniu de măsurare:	2 - 210 µg/L
Instrument:	DR/3000, DR/2000
Numar metoda:	8143
Numar teste:	100
Nume metoda:	Porphyrin
Platforma:	Pernu#e cu pudră