



Silver nitrate digital titrator cartridge, 0.2256 N, 1 piece

Nr. produs:

1439601

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 100 mg/L Cl⁻.

Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 100 mg/L Cl⁻.

For chloride determination by Digital Titrator argentometric titration. Hach method 8207. Digits x 0.1 = mg/L Cl⁻ in 100 mL sample.
Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	10 - 100 mg/L Cl ⁻
Numar metoda:	8207
Numar teste:	40 - 100
Nume metoda:	Silver Nitrate
Parametru:	Chloride
Platforma:	Titrator digital
Titrant:	0.2256 N AgNO ₃