



Sodium Hydroxide Digital Titrator Cartridge, 0.3636 N

Nr. produs:

1437801

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Ready-to-use titrant solution for the determination of Carbon Dioxide, 10-100mg/l CO₂.

Ready-to-use titrant solution for the determination of Carbon Dioxide, 10-100mg/l CO₂.

For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	10 - 100 mg/L CO ₂
Numar metoda:	8205
Numar teste:	40 - 100
Nume metoda:	Phenolphthalein
Parametru:	Carbon Dioxide
Platforma:	Titratore digital
Titrant:	0.3636 N NaOH