




EDTA solution, 1M, 50mL SCDB

Nr. produs:

2241926

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



For determination of calcium and magnesium hardness in water by the Calmagite Colorimetric method, Hach Method 8030. Also used for the determination of calcium and magnesium in ultrapure saturated NaCl feed brine (HNB method). 50 mL self-contained Dropping Bottle.

Specificații

Cantitate:	50 mL
Numar metoda:	8030
Nume metoda:	Calmagite Colorimetric
Parametru:	Calcium and Magnesium Hardness
Platforma:	Solu#ie

Accesorii necesare

- Mănuși de protecție mărimea L, albastre, nitril, fără pudră, 50 de bucăți (Item SM743L)