



Acid Reagent Powder Pillows (pk/100)

Nr. produs:

104299

RON Preț (fără TVA):

Contact

Disponibil

Produse periculoase



Powder reagent in sealed packages for the determination of high range Silica, Silicomolybdate Method. Also used for the determination of high range Hypchlorite, Iodometric Method.

Powder reagent in sealed packages for the determination of high range Silica, Silicomolybdate Method. Also used for the determination of high range Hypchlorite, Iodometric Method.

Specificații

Cantitate:	100 /pk
Numar metoda:	8185
Numar teste:	100
Nume metoda:	Silicomolybdate
Parametru:	Silica
Platforma:	Pernu#e cu pudră