



## Silica Reagent Set, 1-100 mg/L SiO<sub>2</sub>

**Nr. produs:**

**2244300**

RON Preț (fără TVA):  
Disponibil

Contact

Produse periculoase



Powder reagents in sealed packages for the determination of Silica, Silicomolybdate Method. For 25 mL sample. 100 tests.

### Specificații

|                          |                                 |
|--------------------------|---------------------------------|
| Caracteristici speciale: | HR, 25 mL                       |
| Domeniu de măsurare:     | 1 - 100 mg/L SiO <sub>2</sub>   |
| Domeniu de masurare :    | 25 mL                           |
| Instrument:              | DR/3000, DR/2000,<br><br>DR/890 |
| Numar metoda:            | 8185                            |
| Numar teste:             | 100                             |
| Nume metoda:             | Silicomolybdate                 |
| Parametru:               | Silica                          |
| Platforma:               | Pernu#e cu pudră                |