



Mercuric nitrate cartridge for digital titrator, 2.256N

Nr. produs:

92101

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



For titration of chloride, 100 - 8000 mg/L Cl⁻

Specificații

Domeniu de măsurare:	100 - 8000 mg/L Cl ⁻
Numar metoda:	8206
Numar teste:	40 - 100
Nume metoda:	Mercuric Nitrate
Parametru:	Chloride
Platforma:	Titratore digital
Titrant:	2.256 N Hg(NO ₃) ₂