



Chloride Reagent Set, 10-10000mg/L Cl

Nr. produs:

2288000

RON Preț (fără TVA):

Contact

Nu este afișată data expedierii

Produse periculoase



For the determination of chloride by Silver Nitrate Digital Titrator Method (Hach Method 8207). Set contains Chloride 2 Indicator Powder Pillows and two Silver Nitrate Titration Cartridges.

Specificații

Domeniu de măsurare:	10 - 10000 mg/L Cl ⁻
Numar metoda:	8207
Numar teste:	100
Nume metoda:	Silver Nitrate Titration
Platforma:	Titratore digital