



## Reagent powder pillows, SulfaVer 4, sulfate , 25 mL, pk/1000

Nr. produs:

1206528

Produse periculoase

**ARTICOL ÎNVECHIT**

Acest articol nu mai este disponibil.



Produse înlocuitoare sugerate:

- 1206599

Powder reagent in sealed packages for the determination of Sulfate. For 25ml sample.

Powder reagent in sealed packages for the determination of Sulfate. For 25ml sample.

For sulfate determination by the SulfaVer® 4 Turbidimetric Method. USEPA accepted for reporting wastewater analysis using Hach Method 8051. Range: up to 70 mg/L. Sample Size: 25 mL. Pack of 1000 powder pillows.

---

### Specificații

|                      |                                       |
|----------------------|---------------------------------------|
| Cantitate:           | 1000 /pk                              |
| Domeniu de măsurare: | 2 - 70 mg/L SO <sub>4</sub>           |
| EPA compliant:       | Da                                    |
| Instrument:          | DR/4000, DR/3000,<br>DR/2010, DR/2000 |
| Numar metoda:        | 8051                                  |
| Numar teste:         | 1000                                  |
| Nume metoda:         | SulfaVer 4                            |
| Platforma:           | Pernu#e cu pudră                      |