



Reagent set, sulfite, digital titrator

Nr. produs:

2272300

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



For the determination of Sulfite, 4->400mg/l SO₃.

For the determination of Sulfite, 4->400mg/l SO₃.

For determination of sulfite by iodate-iodide Digital Titrator titration, Hach Method 8216. Range: 4 to 400+ mg/L SO₃. Reagent Set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Iodate-Iodide Cartridge (0.3998 N), and Starch Indicator. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Caracteristici speciale:	Sulfite Reagent set, Digital Titrator
Numar metoda:	8216
Numar teste:	100
Nume metoda:	Iodate-Iodide
Platforma:	Titrator digital