



Stablcal Turbidity Standard, 800 NTU, 100 mL

Nr. produs:

2660542

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Ready to use, stabilised Formazine Primary standard solution in 100 mL bottle, 800 NTU.

Ready to use, stabilised Formazine Primary standard solution in 100 mL bottle, 800 NTU.

Stablcal[®] Turbidity Standard for all Turbidimeters. USEPA accepted for calibration purposes. Stablcal[®] Stabilized Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. NTU Value: 800. Bottle: 100 mL.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Specificații

Cantitate:	100 mL
Concentrație:	800 NTU
Platforma:	Solu#ie