



## ULR verification kit, Stablcal, 100 mL

Nr. produs:

2714600

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Standards kit for checking low range turbidities, includes 100 mL of each 0.10 NTU, 0.30 NTU, 0.50 NTU, 1.0 NTU standards, silicone oil and oiling cloth, and instructions.

Standards kit for checking low range turbidities, includes 100 mL of each 0.10 NTU, 0.30 NTU, 0.50 NTU, 1.0 NTU standards, silicone oil and oiling cloth, and instructions.

Stablcal Verification Kit includes: 0.10 NTU, 0.30 NTU, 0.50 NTU, 1.0 NTU standards (100 mL each), silicone oil and oiling cloth, and instructions. For verification only. Not for calibration.

### Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

### Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

### Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

### Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

---

## Specificații

Cantitate: Set of 4

Concentrație: 0,1, 0,3, 0,5, 1,0 NTU

Nota de subsol: \*For verification only. Not for calibration.