



SPADNS 2 Fluoride Reagent, arsenic-free, ACCUVAC, 0.02-2.00 mg/L F

Nr. produs:

2527025

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



ACCUVAC ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Fluoride by the SPADNS 2 method. pk/25

Specificații

Cantitate:	25 /pk
Domeniu de măsurare:	0.02 - 2.00 mg/L F ⁻
EPA compliant:	Da
Method Number 1 Chemistry:	10225
Numar teste:	25
Nume metoda:	SPADNS 2
Parametru:	Fluoride
Platforma:	Accuvac