



Chlorine Standard, 100-130 mg / L 20 mL Amp

Nr. produs:	2835801
RON Preț (fără TVA):	Contact
Disponibil	

This chlorine standard can be diluted down to verify the accuracy of the CL17 chlorine analyser or create a separate calibration curve. Under normal circumstances, a separate calibration curve is not required for the CL17. Package of one 20 mL ampule.

Specificații

Condiții de depozitare:	Keep refrigerated.
Domeniu de măsurare:	N/A
Instrument:	CL17

Accesorii necesare

- Mănuși de protecție mărimea L, albastre, nitril, fără pudră, 50 de bucăți (Item SM743L)