



Stablcal, Turbidity Standard, 40.0 NTU, Bottle/1000 mL

Nr. produs:

2746353

RON Preț (fără TVA):

Contact

Disponibil

Produse periculoase



Stabilised Formazin turbidity Standard for use in any turbidimeter.

Precise, Pre-mixed Formazin Solutions

Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations as low as 0.060 NTU for conventional turbidimeters. Shelf life is guaranteed for two years from the date of manufacture. Whether you are monitoring wastewater effluent, high-quality treated drinking water or an industrial process stream, you can be completely confident in the results you obtain with Hach Stablcal Stabilised Formazin Standards.

No special glassware or difficult lab techniques

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Since then, Formazin has been the only primary standard that is accepted virtually throughout the world for regulatory reporting. All other turbidity standards, both alternate and secondary, must be traced to it.

Specificații

Concentratie: 40 NTU

Condiții de depozitare: Store between 5 - 25 °C, away from direct sunlight

Domeniu de măsurare: N/A

Parametru: Turbidity

Accesorii necesare

- Mănuși de protecție mărimea L, albastre, nitril, fără pudră, 50 de bucăți (Item SM743L)