



Molybdate Reagent Powder Pillows

Nr. produs: 104199

RON Preț (fără TVA): Contact
Disponibil

Powder reagent in sealed package for use in the Silicomolybdate Method for the determination of high range Silica (0-100.0mg/L SiO₂). For 25mL sample. Pack of 100.

Specificații

Cantitate:	100 /pk
Numar metoda:	8185
Numar teste:	100
Nume metoda:	Silicomolybdate
Parametru:	Silica
Platforma:	Pernu#e cu pudră