



Water, deionised and demineralised, 100mL

Nr. produs: 27242
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specificații

Platforma:	Solu#ie
------------	---------