



Water, deionised and demineralised, 500mL

Nr. produs: 27249
RON Preț (fără TVA): Contact
Disponibil

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specificații

Descriere:	Water, Deionized
Platforma:	Solu#ie
Volum:	500 mL