



EDTA Digital Titrator Cartridge, 0.1428 M

Nr. produs:

1496001

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Ready-to-use titrant solution for the analysis of Total Hardness, 1-16 G.d.H. CaCO_3 .

Ready-to-use titrant solution for the analysis of Total Hardness, 1-16 G.d.H. CaCO_3 .

For the determination of Calcium and Total Hardness. 0.1428 M. Digits x 0.1785 = mg/L CaCO_3 in 100 mL, Digits x 0.01 = Gdh.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	1 - 16 G.d.h.
Numar teste:	40 - 100
Nume metoda:	EDTA (ManVer® 2)
Parametru:	Hardness, Total
Platforma:	Titrator digital
Titrant:	0.1428 M EDTA