



Reagent set, hardness (calcium), HR, digital titrator

Nr. produs:

2447500

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



For the determination of Calcium Hardness, 100-4000 mg/L CaCO₃, EDTA Method.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Caracteristici speciale:	Reagent set, HR, Digital Titrator
Domeniu de măsurare:	100 - 4000 mg/L CaCO ₃
Numar metoda:	8204
Numar teste:	40-100
Nume metoda:	EDTA Titration
Parametru:	Hardness (Calcium)
Platforma:	Titrator digital