



## Reagent set, hardness (total), HR, digital titrator

**Nr. produs:**

**2448100**

RON Preț (fără TVA):  
Disponibil

Contact

Produse periculoase



For the determination of Total Hardness, 100-4000 mg/l CaCO<sub>3</sub>, EDTA Method.

For the determination of Total Hardness, 100-4000 mg/l CaCO<sub>3</sub>, EDTA Method.

Reagent set for the determination of Total Hardness by EDTA Digital Titrator titration (Hach method 8213). Range: 100-4,000 mg/L as CaCO<sub>3</sub>. Set contains: Buffer Solution (Hardness 1), ManVer® 2 Hardness Indicator Powder Pillows and 0.800 M EDTA Digital Titrator Cartridge for 40 - 100 tests.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

## Specificații

|                          |                                   |
|--------------------------|-----------------------------------|
| Caracteristici speciale: | Reagent set, HR, Digital Titrator |
| Domeniu de măsurare:     | 100 - 4000 mg/L CaCO <sub>3</sub> |
| Numar metoda:            | 8213                              |
| Numar teste:             | 100                               |
| Nume metoda:             | EDTA Titration                    |
| Parametru:               | Hardness (Total)                  |
| Platforma:               | Titrator digital                  |