



Soluție etalon de clorură de sodiu, 10.000 $\mu\text{S}/\text{cm}$, 500 ml

Nr. produs: 2972249

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 10000 $\mu\text{S}/\text{cm}$ ± 100 $\mu\text{S}/\text{cm}$ at 25°C, 5.556 g/L ± 0.03 g/L as Sodium Chloride (NaCl), (1 $\mu\text{S}/\text{cm}$ \sim 0.5 mg/L TDS), 500 mL bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specificații

Cantitate:	500 mL
Concentrație:	10000 $\mu\text{S}/\text{cm}$ at 25°C ± 100 $\mu\text{S}/\text{cm}$ at 25°C
Condiții de depozitare:	10 - 25 °C (protect from light)
Descriere:	Clorură de sodiu, 10.000 $\mu\text{S}/\text{cm}$, 500 ml
Package settings:	Bottle
Parametru:	Conductivitate
Platforma:	Solu#ie