



Filter paper, glass fiber, extra thick, pore size 1 µm, diameter 47 mm, 100/pk

Nr. produs: 2842100

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

For filtration of heavily contaminated samples. High particulate-holding capacity makes discs efficient depth filters and allows for filtration of large volumes of solutions. Pall Life Sciences.

Specificații

Cantitate:	100
Diametru:	47 mm mm
Dimensiune por:	1 µm
Material:	Glass Microfiber