



## Filter membrane, EZ-Pak, 0.45 micron, 47 mm, black, gridded, 600/pk (millipore)

**Nr. produs:** 2915900

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

---

### Specificații

Cantitate:	600
Diametru:	47 mm mm
Dimensiune por:	0.45 μm
Material:	Mixed-Cellulose Ester