



pH electrode, diameter 12 mm, 10 m fixed moulded cable, pHRET

Nr. produs: LZX476

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Electrode for the determination of the pH value in waste water and pressurised pipes. Optionally with integrated temperature sensor. Suitable for use with the 12 mm inline armature.

Specificații

Adâncime imersie:	205 mm
Condiții de operare:	-5 - 135 °C
Diafragma:	Porous PTFE, double ring diaphragm
Domeniu de măsurare:	0 - 14 pH
Electrolit:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garanție:	6 luni
Interval de presiune:	15 bar
Lungime:	225 mm
Lungime cablu:	10 m fixed cable
Material:	Stainless steel
Reference system:	Ag/AgCl
Senzor de temperatura:	None