



pH glass electrode, refillable, pure water applications

Nr. produs:	LZY025
RON Preț (fără TVA):	Contact
Disponibil	

pH electrode for use in ultra-pure water

Specially developed for use in ultra-pure water

Refillable glass electrodes for special requirements and high-precision measurements

Improved durability through triple ceramic diaphragm

Easy to install in stainless steel or polypropylene flow fittings

Specificații

Condiții de operare:	-30 - 80 °C
Conexiunea cablului:	Plug head S7
Diafragma:	Triple ceramic
Domeniu de măsurare:	0 - 12 pH
Electrolit:	KCl liquid
Garanție:	6 luni
Greutate:	0,15 kg
Lungime:	120 mm
Reference system:	Ag/AgCl
Senzor de temperatura:	None
Timp de raspuns T98%:	< 30 s
Tip electrod:	Glass