



Pole mounting hardware pH, 24cm bracket, PVC pole 2m

Nr. produs: LZY714.99.22910

RON Preț (fără TVA): Contact

Expediere în 1 săptămână

PVC pole mounting hardware for 1" NPT pH sensor including 24 cm SS bracket and 2 m PVC pole. Suitable for all 1" NPT pH sensors like 1200S sc (LXV426.99.10001) and pHD ORP (LXV427.99.10001)

Specificații

Diametru:	40 mm
Greutate:	Lifting part only: 2.5 kg
Instrument:	Suitable for all 1" NPT pH sensors
Lungime:	Pole: 2000 mm
Material:	Wetted parts:
	pole PVC;
	sensor adapter PVC;
	screw set SS.